## **Purism FFP2 Face Masks**



GB 2626-2006 CE 729899 + 729900 BS EN 149:2001 +A1:2009 EU 2016/425 EU 2020/403 FDA GTS22151

# GUARANTEED GERTIFIED Independently Tested & Verified

















#### **Purism FFP2 Face Masks**





## **Purism FFP2 Face Masks**





#### Brand : Purism

Product Name : KN95 FFP2 Protective Face Mask

**Product Performance** : Perfect tight fit for all face shapes providing optimal protection. Suitable for protection against non-oily particles such as bacteria, flu virus, dust, pollen etc

**Fibre Fabric Construction** : Spun-bonded, nonwoven fabric and melt-blown non-woven fabric with sonic welding

**Application** : Schools, airplanes, public transport, workplaces and public areas where social distancing cannot be achieved.

## **Five Layers of Filtration**

Two layers of spun-bonded, non-woven fabric, two layers of melt-blown, non-woven fabric, one layer of air-blown, non-woven fabric, to effectively filter out viruses and bacteria to protect your health and the health of those around you.



- Spun-bond, non-woven
- 3 Air-blown, non-woven
- Spun-bon, non-woven
- S Extra-soft earloops

- Melt-blown, non-woven
- 4 Melt-blown, non-woven
- 6 Three-dimensional nose-bridge strip

**Testing bodies** 

**Certificate date** 

Classification Certificate number



#### HSE safety warn against fake KN95 face masks

There have been a 'substantial' number of face masks claiming to be of KN95 standards available in the UK over the course of the COVID crisis which provide an **inadequate level of protection**.



The Health and Safety Executive have stated there are many masks of poor quality which are accompanied by fake or fraudulent paperwork.

When purchasing face masks and filtering facepiece (FFP) respirators you should be assured that these are fit for purpose and your supplier should be able to provide certification which can be easily verified with the governing and/or 'Notified' testing body.

Masks or mask packaging should be CE marked and the CE certification should be verifiable.

FFP2 masks should be tested to the European standard BSEN149:2001+A1:2009

Mask testing should be accompanied by Notified Body documentation showing compliance with essential health and safety requirements as required by the Personal Protective Equipment Regulations (EU) 2016/425.